
**BiomarC is currently available through a distribution
agreement with Vigeo Healthcare Solutions**

Contact your Vigeo representative to order today.

Tel. +39 0376 449097

Fax +39 0376 282245

Catalog Number	Compatible With	Shape, Size	(Modality) Applicator Description
040300	Core Biopsy Needles	Petite Barbell, 1x3mm	(U/S) Rigid SS 14g needle
040320		Tribell, 1x5mm	(U/S) Rigid SS 14g needle
040306		Barbell, 2x4mm	(U/S) Rigid SS 13g needle
040310		Tribell, 1x5mm	(U/S) Rigid SS 17g needle
040204	Mammotome elite 10g (through introducer)	Barbell, 2x4mm	(U/S) flexible blunt end-deploy applicator
040203	Mammotome 11g	Barbell, 2x4mm	(ST, U/S) flexible, side-deploy applicator
040220	Mammotome 8g	Barbell, 2x4mm	(ST, U/S) flexible, side-deploy applicator
040204	Hologic Eviva 9g	Barbell, 2x4mm	(ST) flexible, blunt end-deploy applicator
040203	Hologic ATEC 9g (detached needle)	Barbell, 2x4mm	(ST, U/S) flexible, side-deploy applicator
040228	Hologic ATEC 9g (back of probe)	Barbell, 2x4mm	(ST, U/S) flexible, side-deploy applicator
040239	Bard ENCOR 10g	Barbell, 2x4mm	(ST, U/S) flexible, side-deploy applicator



Carbon Medical Technologies, Inc
1290 Hammond Road
Saint Paul, Minnesota 55110 USA
Phone: 888-207-0262
Fax: 651-407-1975
e-mail: info@carbonmed.com
www.carbonmed.com

©2021 Carbon Medical Technologies, Inc.
All rights reserved.
BiomarC is a registered trademark of
Carbon Medical Technologies, Inc.